

Product Number: **RDS0004**
 Product Description: Opiates L1-L3 (Multilevel)
URINE TOXICOLOGY CONTROL

Lot Number: C8914		Expiration Date: 2023-06-30		
Analyte	Target Value	Verified Value	Units	Thawed Stability Days
Codeine	7500	6000	ng/mL	25
Dihydrocodeine	7500	6400	ng/mL	25
Hydrocodone	7500	5200	ng/mL	25
Norhydrocodone	7500	5331	ng/mL	25
Hydromorphone	7500	5500	ng/mL	25
Morphine	7500	6600	ng/mL	25
Naloxone	7500	6300	ng/mL	25
Oxycodone	7500	4800	ng/mL	25
Noroxycodone	7500	5025	ng/mL	25
Oxymorphone	7500	6100	ng/mL	25
Noroxymorphone	7500	5791	ng/mL	25

Lot Number: C8913		Expiration Date: 2023-06-30		
Analyte	Target Value	Verified Value	Units	Thawed Stability Days
Codeine	2000	1800	ng/mL	25
Dihydrocodeine	2000	1900	ng/mL	25
Hydrocodone	2000	1500	ng/mL	25
Norhydrocodone	2000	1638	ng/mL	25
Hydromorphone	2000	1600	ng/mL	25
Morphine	2000	1700	ng/mL	25
Naloxone	2000	1700	ng/mL	25
Oxycodone	2000	1600	ng/mL	25
Noroxycodone	2000	1358	ng/mL	25
Oxymorphone	2000	1600	ng/mL	25
Noroxymorphone	2000	1708	ng/mL	25

Lot Number: C8912		Expiration Date: 2023-06-30		
Analyte	Target Value	Verified Value	Units	Thawed Stability Days
Codeine	35	30	ng/mL	25
Dihydrocodeine	35	33	ng/mL	25
Hydrocodone	35	28	ng/mL	25
Norhydrocodone	35	N/A	ng/mL	25
Hydromorphone	35	28	ng/mL	25
Morphine	35	35	ng/mL	25
Naloxone	35	30	ng/mL	25
Oxycodone	35	31	ng/mL	25
Noroxycodone	35	N/A	ng/mL	25
Oxymorphone	35	31	ng/mL	25
Noroxymorphone	35	N/A	ng/mL	25

Product Number: RDS0004
Product Description: Opiates L1-L3 (Multilevel)
URINE TOXICOLOGY CONTROL

PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For in vitro diagnostic use only.
3. For analytical use only.

PROCEDURE:

1. Allow fresh frozen control material to thaw at room temperature with cap on.
2. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
3. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.
4. Assay control material in same manner as patient specimens, following the exact same instructions from the entire test method.
5. Record the results obtained on a quality control chart that describes statistical limits for the test method and the particular lot of control material.

STORAGE AND STABILITY:

1. Store fresh frozen control material at or below -10°C (14°F). Stable to expiration date printed on the insert and label.
2. Store thawed control material at or below 2-8°C (35-46°F). **See table for thawed stability.**

For technical assistance email UTAK's R&D Department at <awards@utak.com>

UTAK Laboratories, Inc.

25020 Avenue Tibbitts
Valencia, CA 91355
USA

T: (888) 882-5522

F: (661) 294-9272

Product Number: RDS0004
Product Description: Opiates L1-L3 (Multilevel)
Lot Number: C8912-C8914
Expiration Date: 2023-06-30
Volume: 1x1mL Vials, Frozen