

PRECAUTIONS:

1. Although the serum has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For in research use only.

STORAGE AND STABILITY:

1. Store fresh frozen control material at or below -10°C (14°F). Stable to expiration date printed on the insert and label.
2. Store thawed control material at or below 2-8°C (35-46°F).

PROCEDURE:

1. Allow fresh frozen control material to thaw at room temperature with cap on.
2. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
3. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.
4. Assay control material in same manner as patient specimens, following the exact same instructions from the entire test method.
5. Record the results obtained on a quality control chart that describes statistical limits for the test method and the particular lot of control material.

EXPECTED VALUES:

1. Laboratories should establish their own statistical values for precision and expected ranges.

Lot Number:	RD255	Expiration Date:	2024-11-30
	Analyte	Target Value	Units
	Fluconazole	25000	ng/mL
	Hydroxy-Itraconazole	2500	ng/mL
	Aripiprazole	2500	ng/mL
	Posaconazole	5000	ng/mL
	Voriconazole	7500	ng/mL

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

For technical assistance contact UTAK's R&D Department at: awards@utak.com

UTAK Laboratories, Inc.
 25020 Avenue Tibbitts
 Valencia, CA 91355
 USA
 T: (888) 882-5522
 F: (661) 294-9272
 E: welovecontrol@utak.com

Product Number: RDS0005
Product Description: Antifungals
Lot Number: RD255
Expiration Date: 2024-11-30
Volume: 1x1mL Vials, Frozen

